

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10558 - IMMUNIZING PEPTIDE

SKU: EBP10558

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

ID

NP 005268.1; NP 963885.1; NP 963886.1

Peptide with sequence C-QTYFEMARTAGD, from the internal region of the protein sequence according to Immunogen

NP 005268.1; NP 963885.1; NP 963886.1.

The peptide was used in the production of EB10558 - Goat Anti-GPM6A Antibody. It is currently untested in **Product**

Comments blocking applications. **Peptide**

C-QTYFEMARTAGD Sequence

Shipping

Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png