

EBP10581 - IMMUNIZING PEPTIDE

SKU: EBP10581



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_000388.2
Immunogen	Peptide with sequence CDKQEFQDAWKLPK, from the internal region of the protein sequence according to NP_000388.2.
Product	The peptide was used in the production of EBP10581 - Goat Anti-GP91-PHOX / NOX2 Antibody. It is currently
Comments	untested in blocking applications.
Peptide Sequence	CDKQEFQDAWKLPK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png