

EBP10594 - IMMUNIZING PEPTIDE

SKU: EBP10594



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_005260.1; NP_001153517.1; NP_001161081.1

Immunogen Peptide with sequence TYSQCPRLEHYGQ, from the internal region of the protein sequence according to NP_005260.1; NP_001153517.1; NP_001161081.1.

Product Comments The peptide was used in the production of EBP10594 - Goat Anti-GLI1 Antibody. It is currently untested in blocking applications.

Peptide Sequence TYSQCPRLEHYGQ

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png