

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10601 - IMMUNIZING PEPTIDE

SKU: EBP10601

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_001028300.1 ID

 $\begin{tabular}{ll} \textbf{Immunogen} & Peptide with sequence C-NGSIYILRRDSKHPK, from the internal region of the protein sequence according to NP_001028300.1. \end{tabular}$

The peptide was used in the production of EB10601 - Goat Anti-Bbs1 (mouse) Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide C-NGSIYILRRDSKHPK

Sequence

Ambient Instructions Shipping

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade