

EBP10619 - IMMUNIZING PEPTIDE

SKU: EBP10619



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_997126.2

Immunogen Peptide with sequence ECKRSVKYNYETFQ, from the C Terminus of the protein sequence according to NP_997126.2.

Product Comments The peptide was used in the production of EBP10619 - Goat Anti-Mucin 19 / Apomucin Antibody. It is currently untested in blocking applications.

Peptide Sequence ECKRSVKYNYETFQ

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

DOCUMENTS

- [Data Sheet](#)