



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP10619 - IMMUNIZING PEPTIDE**

**SKU:** EBP10619

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage Shipped at ambient temperature, store at -20°C.

**Accession** 

NP\_997126.2 ID

Immunogen Peptide with sequence ECKRSVKYNYETFQ, from the C Terminus of the protein sequence according to NP\_997126.2.

**Product** The peptide was used in the production of EB10619 - Goat Anti-Mucin 19 / Apomucin Antibody. It is currently

**Comments** untested in blocking applications.

**Peptide** 

**ECKRSVKYNYETFQ** 

Sequence

Shipping Ambient Instructions

**Product** Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite\_plus\_medium.png$