

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10642 - IMMUNIZING PEPTIDE

SKU: EBP10642

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 005620.1; NP 001136277.1; NP 001136278.1 ID

Peptide with sequence ASLANLTCDQLADRR, from the internal region of the protein sequence according to Immunogen

NP_005620.1; NP_001136277.1; NP_001136278.1.

The peptide was used in the production of EB10642 - Goat Anti-Creatine transporter 1 / SLC6A8 Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

ASLANLTCDQLADRR Sequence

Shipping

Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png