

EBP10666 - IMMUNIZING PEPTIDE

SKU: EBP10666



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_006437.3

Immunogen Peptide with sequence C-KYQEKDINASEN, from the internal region (near C Terminus) of the protein sequence according to NP_006437.3.

Product Comments The peptide was used in the production of EBP10666 - Goat Anti-SLCO1B1 Antibody. It is currently untested in blocking applications.

Peptide Sequence C-KYQEKDINASEN

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png