

EBP10670 - IMMUNIZING PEPTIDE

SKU: EBP10670



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_060664.2; NP_001164035.1; NP_001164059.1; NP_001164058.1
Immunogen	Peptide with sequence C-KQMQQQEHDSDKAS, from the internal region of the protein sequence according to NP_060664.2; NP_001164035.1; NP_001164059.1; NP_001164058.1.
Product Comments	The peptide was used in the production of EBP10670 - Goat Anti-HHAT (aa227-240) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-KQMQQQEHDSDKAS
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)