

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10681 - IMMUNIZING PEPTIDE

SKU: EBP10681

SPECIFICATIONS

Formulation $100\mu g$ of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

Accession

NP 001146956.1; NP 001147135.1; NP 001147024.1; NP 115936.2; NP 001147107.1; NP 001147162.1 ID

Peptide with sequence C-HTLEDLSARVDAVKE, from the internal region (near C Terminus) of the protein

Immunogen sequence according to NP 001146956.1; NP 001147135.1; NP 001147024.1; NP 115936.2; NP 001147107.1;

NP 001147162.1.

The peptide was used in the production of EB10681 - Goat Anti-SCOC Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

C-HTLEDLSARVDAVKE Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade