

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10693 - IMMUNIZING PEPTIDE

SKU: EBP10693

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Instructions Shipped at ambient temperature, store at -20°C.

ID

Accession $NP_004405.3;\ NP_981962.1;\ NP_981963.1;\ NP_001272320.2;\ NP_001272318.1;\ NP_001317945.1$

Peptide with sequence C-HIGSSHLAPPNPD, from the internal region of the protein sequence according to

Immunogen NP_004405.3; NP_981962.1; NP_981963.1; NP_001272320.2; NP_001272318.1; NP_001317945.1.

The peptide was used in the production of EB10693 - Goat Anti-Calcipressin-1 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

C-HIGSSHLAPPNPD Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade