

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10696 - IMMUNIZING PEPTIDE

SKU: EBP10696

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 000389.1; NP 001123291.1; NP 001165131.1 ID

Peptide with sequence CLPNLDHVGHPTER, from the C Terminus of the protein sequence according to Immunogen

 $NP_000389.1;\ NP_001123291.1;\ NP_001165131.1.$

The peptide was used in the production of EB10696 - Goat Anti-CYB5R3 / Dia 1 Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

CLPNLDHVGHPTER Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade