

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP10738 - IMMUNIZING PEPTIDE**

**SKU:** EBP10738

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_787058.5; NP\_008820.4; NP\_001193853.2; NP\_001193854.2 ID

**Immunogen** Peptide with sequence KKRPRNSYKCGTNT, from the internal region of the protein sequence according to NP\_787058.5; NP\_008820.4; NP\_001193853.2; NP\_001193854.2.

The peptide was used in the production of EB10738 - Goat Anti-CD86 / B7-2 (aa273-286) Antibody. It is currently Product

**Comments** untested in blocking applications.

Peptide

KKRPRNSYKCGTNT Sequence

Shipping

Instructions Ambient

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade