

EBP10738 - IMMUNIZING PEPTIDE

SKU: EBP10738



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_787058.5; NP_008820.4; NP_001193853.2; NP_001193854.2
Immunogen	Peptide with sequence KKRPRNSYKCGTNT, from the internal region of the protein sequence according to NP_787058.5; NP_008820.4; NP_001193853.2; NP_001193854.2.
Product Comments	The peptide was used in the production of EB10738 - Goat Anti-CD86 / B7-2 (aa273-286) Antibody. It is currently untested in blocking applications.
Peptide Sequence	KKRPRNSYKCGTNT
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png