

EBP10739 - IMMUNIZING PEPTIDE

SKU: EBP10739



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_787058.5; NP_008820.4; NP_795711.2
Immunogen	Peptide with sequence C-EKFDSVHSKYMGRS, from the internal region of the protein sequence according to NP_787058.5; NP_008820.4; NP_795711.2.
Product	The peptide was used in the production of EBP10739 - Goat Anti-CD86 / B7-2 (aa 73-87) Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-EKFDSVHSKYMGRS
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png