

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP10746 - IMMUNIZING PEPTIDE**

**SKU:** EBP10746

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_001028508.2; NP\_001157147.1; NP\_001157148.1 ID

Peptide with sequence RAKKDSRKPSLTK-C, from the N Terminus (near) of the protein sequence according to Immunogen

NP\_001028508.2; NP\_001157147.1; NP\_001157148.1.

The peptide was used in the production of EB10746 - Goat Anti-Abcc4 (mouse) Antibody. It is currently untested in Product

Comments blocking applications.

**Peptide** 

**RAKKDSRKPSLTK-C** Sequence

Instructions Ambient Shipping

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade