

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10772 - IMMUNIZING PEPTIDE

SKU: EBP10772

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

NP_000894.1 ID

 $\begin{tabular}{ll} \textbf{Immunogen} \\ \textbf{NP_000894.1}. \\ \end{tabular} \begin{tabular}{ll} \textbf{Peptide with sequence C-DKGPFRSKKAVLS, from the internal region of the protein sequence according to NP_000894.1}. \\ \end{tabular}$

Product The peptide was used in the production of EB10772 - Goat Anti-NQO1 (isoform a) Antibody. It is currently untested

Comments in blocking applications.

Peptide

C-DKGPFRSKKAVLS Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade