

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10779 - IMMUNIZING PEPTIDE

SKU: EBP10779

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 944494.1; NP 000807.2; NP 944493.2 ID

Peptide with sequence KSDDDYEDYASNK-C, from the N Terminus (mature) of the protein sequence according Immunogen to NP_944494.1; NP_000807.2; NP_944493.2.

The peptide was used in the production of EB10779 - Goat Anti-GABRG2 Antibody. It is currently untested in **Product**

Comments blocking applications. **Peptide**

KSDDDYEDYASNK-C Sequence

Shipping Instructions Ambient

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade