

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10782 - IMMUNIZING PEPTIDE

SKU: EBP10782

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_000546.2 ID

 $\begin{tabular}{l|l} \textbf{Immunogen} & Peptide with sequence DRSHCQSLRFHQN, from the internal region of the protein sequence according to NP_000546.2. \end{tabular}$

Product The peptide was used in the production of EB10782 - Goat Anti-Doublecortin / DCX (aa69-81) Antibody. It is

Comments currently untested in blocking applications.

Peptide

DRSHCQSLRFHQN Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade