

EBP10789 - IMMUNIZING PEPTIDE

SKU: EBP10789



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001521.1; NP_851397.1
Immunogen	Peptide with sequence C-KMTQLFTKVESED, from the internal region of the protein sequence according to NP_001521.1; NP_851397.1.
Product	The peptide was used in the production of EBP10789 - Goat Anti-HIF1alpha / MOP1 Antibody. It is currently
Comments	untested in blocking applications.
Peptide Sequence	C-KMTQLFTKVESED
Shipping	
Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png