

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10837 - IMMUNIZING PEPTIDE

SKU: EBP10837

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_997125.2 ID

The peptide was used in the production of EB10837 - Goat Anti-Npc1l1 (mouse) (aa56-69) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

Sequence

CLSNTPARHVTGDH

Shipping

Ambient Instructions

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png