

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10841 - IMMUNIZING PEPTIDE

SKU: EBP10841

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 001005242.2; NP 004563.2 ID

Peptide with sequence C-RGTAQYSSQKSVEER, from the internal region of the protein sequence according to Immunogen

NP 001005242.2; NP 004563.2.

The peptide was used in the production of EB10841 - Goat Anti-Plakophilin 2 (aa125-139) Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

C-RGTAQYSSQKSVEER Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade