

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10844 - IMMUNIZING PEPTIDE

SKU: EBP10844

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

 $NP_932173.1; NP_000326.2; NP_001092874.1; NP_001092875.1; NP_001153632.1; NP_001153633.1; NP_001155633.1; NP_001155633.1; NP_001155633.1; NP_001155633.1; NP_001155635.1; NP_001155655.1; NP_001155655.1; NP_001155655.1; NP_001155655.1; NP_001155655.1; NP_001155655.1; NP_00155655.1; NP_00155655.1; NP_00155655.1; NP_00155655.1; NP_001$ ID

Immunogen Peptide with sequence C-PTRKETRFEEGEQP, from the internal region of the protein sequence according to NP_932173.1; NP_000326.2; NP_001092874.1; NP_001092875.1; NP_001153632.1; NP_001153633.1.

The peptide was used in the production of EB10844 - Goat Anti-SCN5A (aa1021-1034) Antibody. It is currently Product

Comments untested in blocking applications.

Peptide

C-PTRKETRFEEGEQP Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade