

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10857 - IMMUNIZING PEPTIDE

SKU: EBP10857

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 072046.2; NP 715640.2; NP 783329.1 ID

Peptide with sequence C-DKNATKAGDDEPEY, from the internal region of the protein sequence according to Immunogen

NP 072046.2; NP 715640.2; NP 783329.1.

The peptide was used in the production of EB10857 - Goat Anti-DNMT3A Antibody. It is currently untested in **Product**

Comments blocking applications. **Peptide**

Sequence

C-DKNATKAGDDEPEY

Shipping

Instructions Ambient

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png