

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10862 - IMMUNIZING PEPTIDE

SKU: EBP10862

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

100 μg **Unit Size**

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 006068.2; NP 001182447.1 ID

Peptide with sequence C-EYENRIRAYSTPDK, from the internal region of the protein sequence according to Immunogen

NP 006068.2; NP 001182447.1.

The peptide was used in the production of EB10862 - Goat Anti-CBARA1 (aa113-126) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-EYENRIRAYSTPDK Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade