

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10872 - IMMUNIZING PEPTIDE

SKU: EBP10872

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_006151.3 ID

Peptide with sequence C-DGEDDEPRSDKGDAP, from the internal region (near N Terminus) of the protein sequence

according to NP_006151.3.

The peptide was used in the production of EB10872 - Goat Anti-NEUROD2 (aa23-37) Antibody. It is currently Product

Comments untested in blocking applications.

Peptide

C-DGEDDEPRSDKGDAP Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade