

EBP10879 - IMMUNIZING PEPTIDE

SKU: EBP10879



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_110437.2; NP_001139021.1
Immunogen	Peptide with sequence C-RDLESLREYVESQLQ, from the internal region of the protein sequence according to NP_110437.2; NP_001139021.1.
Product Comments	The peptide was used in the production of EBP10879 - Goat Anti-TXNDC5 (aa281-295) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-RDLESLREYVESQLQ
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)