

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP10880-B - IMMUNIZING PEPTIDE**

**SKU:** EBP10880-B

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

100 μg **Unit Size** 

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 006840.2 ID

Peptide with sequence C-EYKSTRDLETFSK., from the internal region (near C Terminus) of the protein sequence **Immunogen** 

according to NP 006840.2.

The peptide was used in the production of EB10880-B - Goat Anti-PDIA2 (aa477-481), Biotinylated Antibody. It **Product** 

**Comments** is currently untested in blocking applications.

**Peptide** 

C-EYKSTRDLETFSK. Sequence

**Shipping** 

Instructions Ambient

**Product** https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade