

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10888 - IMMUNIZING PEPTIDE

SKU: EBP10888

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

 $NP_000255.2; \ NP_001077072.1; \ NP_001077071.1; \ NP_001077073.1$ ID

Peptide with sequence C-EELKKINYQVDSWEE, from the internal region of the protein sequence according to Immunogen NP_000255.2; NP_001077072.1; NP_001077071.1; NP_001077073.1.

The peptide was used in the production of EB10888 - Goat Anti-PTCH1 (aa266-280) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-EELKKINYQVDSWEE Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade