

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10898 - IMMUNIZING PEPTIDE

SKU: EBP10898

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession

 $NP_038474.2;\ NP_001138597.1;\ NP_001278592.1$ ID

 $\textbf{Immunogen} \begin{tabular}{l} \textbf{Peptide with sequence C-RYEHSKPLKQEEAT, from the internet regoin of the protein sequence according to NP_038474.2; NP_001138597.1; NP_001278592.1. \end{tabular}$

The peptide was used in the production of EB10898 - Goat Anti-MKRN1 (aa105-118) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-RYEHSKPLKQEEAT Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$ Grade