

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10922 - IMMUNIZING PEPTIDE

SKU: EBP10922

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

100 μg **Unit Size**

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 061216.2; NP 001152771.1 ID

Peptide with sequence C-KPEGQEQFQVAEKP, from the internal region of the protein sequence according to Immunogen

NP 061216.2; NP 001152771.1.

The peptide was used in the production of EB10922 - Goat Anti-Itih4 (mouse) Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

C-KPEGQEQFQVAEKP Sequence

Shipping

Instructions Ambient

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png