

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10930 - IMMUNIZING PEPTIDE

SKU: EBP10930

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

100 μg **Unit Size**

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 038754.1; NP 001177646.1 ID

Peptide with sequence C-DKKRTEYLPAHEDRT, from the internal region of the protein sequence according to Immunogen

NP_038754.1; NP_001177646.1.

The peptide was used in the production of EB10930 - Goat Anti-Dbf4 (mouse) Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide C-DKKRTEYLPAHEDRT

Sequence

Shipping Instructions Ambient

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade