

EBP10933 - IMMUNIZING PEPTIDE

SKU: EBP10933



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_056998.4
Immunogen	Peptide with sequence C-QLDYAEYVRRDSE, from the internal region of the protein sequence according to NP_056998.4.
Product Comments	The peptide was used in the production of EBP10933 - Goat Anti-TXNDC11 (aa113-126) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-QLDYAEYVRRDSE
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)