

EBP10945 - IMMUNIZING PEPTIDE

SKU: EBP10945



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_112473.1
Immunogen	Peptide with sequence KQAYEEPRQDHELK-C, from the N Terminus of the protein sequence according to NP_112473.1.
Product Comments	The peptide was used in the production of EBP10945 - Goat Anti-SLC19A1 / RFC1 (mouse) Antibody. It is currently untested in blocking applications.
Peptide Sequence	KQAYEEPRQDHELK-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)