

EBP10954-B - IMMUNIZING PEPTIDE

SKU: EBP10954-B



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_055121.1
Immunogen	Peptide with sequence C-QADPNKVSARAKKR., from the internal region of the protein sequence according to NP_055121.1.
Product	The peptide was used in the production of EBP10954-B - Goat Anti-HSPC117 (aa201-215), Biotinylated Antibody.
Comments	It is currently untested in blocking applications.
Peptide Sequence	C-QADPNKVSARAKKR.
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)