

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10954 - IMMUNIZING PEPTIDE

SKU: EBP10954

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

ID

NP_055121.1

Immunogen Peptide with sequence C-QADPNKVSARAKKR, from the internal region of the protein sequence according to NP_055121.1.

The peptide was used in the production of EB10954 - Goat Anti-HSPC117 (aa201-215) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

Sequence

C-QADPNKVSARAKKR

Shipping

Ambient Instructions

Product

Grade

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$