

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10972 - IMMUNIZING PEPTIDE

SKU: EBP10972

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 000627.2; NP 001019637.1 ID

Peptide with sequence C-EAIKRDFGSFDK, from the internal region of the protein sequence according to Immunogen

NP 000627.2; NP 001019637.1.

The peptide was used in the production of EB10972 - Goat Anti-MNSOD (aa119-130) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide C-EAIKRDFGSFDK

Sequence

Shipping Instructions Ambient

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade