

# EBP11005 - IMMUNIZING PEPTIDE

**SKU:** EBP11005



---

## SPECIFICATIONS

|                              |                                                                                                                                                                                                                                         |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Formulation</b>           | 100µg of Lyophilised peptide, formulation De-ionized (DI) water.                                                                                                                                                                        |
| <b>Unit Size</b>             | 100 µg                                                                                                                                                                                                                                  |
| <b>Storage Instructions</b>  | Shipped at ambient temperature, store at -20°C.                                                                                                                                                                                         |
| <b>Accession ID</b>          | NP_000515.2                                                                                                                                                                                                                             |
| <b>Immunogen</b>             | Peptide with sequence C-QPKKSVNGESGRN, from the internal region of the protein sequence according to NP_000515.2.                                                                                                                       |
| <b>Product Comments</b>      | The peptide was used in the production of EB11005 - Goat Anti-HTR1A (aa249-262) Antibody. It is currently untested in blocking applications.                                                                                            |
| <b>Peptide Sequence</b>      | C-QPKKSVNGESGRN                                                                                                                                                                                                                         |
| <b>Shipping Instructions</b> | Ambient                                                                                                                                                                                                                                 |
| <b>Product Grade</b>         | <a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png</a> |

## DOCUMENTS

- [Data Sheet](#)