

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP11019 - IMMUNIZING PEPTIDE**

**SKU:** EBP11019

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_001010973.1; NP\_032312.2 ID

 $\label{eq:lower_lower_lower} \textbf{Immunogen} \ \ \frac{\text{Peptide with sequence C-NRNSHKTSLRSNASQ, from the internal region of the protein sequence according to NP_001010973.1; NP_032312.2.}$ 

The peptide was used in the production of EB11019 - Goat Anti-Histamine Receptor H2 (mouse) Antibody. It is Product

**Comments** currently untested in blocking applications.

**Peptide** 

C-NRNSHKTSLRSNASQ Sequence

Shipping

Ambient Instructions

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade