

EBP11042 - IMMUNIZING PEPTIDE

SKU: EBP11042



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_006155.2; NP_001138884.1; NP_001138885.1; NP_001300831.1;
Immunogen	Peptide with sequence C-KHSSRLEAHLTRDE, from the internal region of the protein sequence according to NP_006155.2; NP_001138884.1; NP_001138885.1; NP_001300831.1;.
Product Comments	The peptide was used in the production of EB11042 - Goat Anti-NRF2 (aa445-458) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-KHSSRLEAHLTRDE
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

DOCUMENTS

- [Data Sheet](#)