



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11042 - IMMUNIZING PEPTIDE

SKU: EBP11042

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

NP_006155.2; NP_001138884.1; NP_001138885.1; NP_001300831.1; ID

Peptide with sequence C-KHSSRLEAHLTRDE, from the internal region of the protein sequence according to Immunogen NP_006155.2; NP_001138884.1; NP_001138885.1; NP_001300831.1;.

The peptide was used in the production of EB11042 - Goat Anti-NRF2 (aa445-458) Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

C-KHSSRLEAHLTRDE Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade