

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11060 - IMMUNIZING PEPTIDE

SKU: EBP11060

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_031426.2 ID

Immunogen Peptide with sequence KKTTMEKLRCVHP, from the internal region of the protein sequence according to NP_031426.2.

The peptide was used in the production of EB11060 - Goat Anti-Adam12 (mouse) Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

KKTTMEKLRCVHP Sequence

Ambient Instructions Shipping

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade