

EBP11063 - IMMUNIZING PEPTIDE

SKU: EBP11063



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_000603.1; NP_001121070.1

Immunogen Peptide with sequence ETYCATPAKSERD, from the internal region of the protein sequence according to NP_000603.1; NP_001121070.1.

Product Comments The peptide was used in the production of EBP11063 - Goat Anti-IGF2 (aa81-93) Antibody. It is currently untested in blocking applications.

Peptide Sequence ETYCATPAKSERD

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)