

EBP11070 - IMMUNIZING PEPTIDE

SKU: EBP11070



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_061968.1; NP_001120930.1
Immunogen	Peptide with sequence C-RSALGSALDPLADK, from the internal region of the protein sequence according to NP_061968.1; NP_001120930.1.
Product Comments	The peptide was used in the production of EBP11070 - Goat Anti-cardiolipin synthase 1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-RSALGSALDPLADK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)