



EBP11073 - IMMUNIZING PEPTIDE

SKU: EBP11073



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_031869.1

Immunogen Peptide with sequence TKLQGTCKPDKPN, from the C Terminus of the protein sequence according to NP_031869.1.

Product The peptide was used in the production of EBP11073 - Goat Anti-defensin, beta 1 (mouse) Antibody. It is currently

Comments untested in blocking applications.

Peptide Sequence TKLQGTCKPDKPN

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)