

EBP11077 - IMMUNIZING PEPTIDE

SKU: EBP11077



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_000761.3; NP_001185782.1; NP_001185783.1
Immunogen	Peptide with sequence C-KNVALTRSYIREK, from the internal region of the protein sequence according to NP_000761.3; NP_001185782.1; NP_001185783.1.
Product Comments	The peptide was used in the production of EBP11077 - Goat Anti-cytochrome P450 2C8 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-KNVALTRSYIREK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)