

EBP11080 - IMMUNIZING PEPTIDE

SKU: EBP11080



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_004604.2; NP_945189.1

Immunogen Peptide with sequence NGRDHHTADLCRE, from the N Terminus of the protein sequence according to NP_004604.2; NP_945189.1.

Product Comments The peptide was used in the production of EBP11080 - Goat Anti-Transglutaminase 2 (aa17-29) Antibody. It is currently untested in blocking applications.

Peptide Sequence NGRDHHTADLCRE

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png