

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP11081 - IMMUNIZING PEPTIDE**

**SKU:** EBP11081

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 004604.2; NP 001310245.1; NP 001310247.1; NP 001310246.1 ID

Peptide with sequence C-ERDLYLENPEIKIR, from the internal region of the protein sequence according to Immunogen

NP 004604.2; NP 001310245.1; NP 001310247.1; NP 001310246.1.

The peptide was used in the production of EB11081 - Goat Anti-Transglutaminase 2 (aa578-591) Antibody. It is **Product** 

**Comments** currently untested in blocking applications.

**Peptide** 

C-ERDLYLENPEIKIR Sequence

**Shipping** 

Instructions Ambient

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade