

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11085 - IMMUNIZING PEPTIDE

SKU: EBP11085

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 003232.2 ID

Peptide with sequence C-NQPLQSYHQLKLEFS, from the internal region (near N Terminus) of the protein Immunogen

sequence according to NP 003232.2.

The peptide was used in the production of EB11085 - Goat Anti-transglutaminase 4 (aa49-63) Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

C-NQPLQSYHQLKLEFS Sequence

Shipping

Instructions Ambient

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade