

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11087 - IMMUNIZING PEPTIDE

SKU: EBP11087

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

100 μg **Unit Size**

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 780681.2 ID

Peptide with sequence C-QKSAEKQGDLTPSET, from the internal region of the protein sequence according to **Immunogen**

NP 780681.2.

The peptide was used in the production of EB11087 - Goat Anti-Terminal Uridylyltransferase 4 (mouse) **Product**

Comments Antibody. It is currently untested in blocking applications.

Peptide

C-QKSAEKQGDLTPSET Sequence

Shipping

Grade

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png